

Pushbutton (BTN) Sensor Data Sheet

BTN 21072016

SPECIFICATIONS

> Type: Momentary

FEATURES

> Ergonomic design
> Easy-to-use

APPLICATIONS

> Event annotation

GENERAL DESCRIPTION

A common need when recording biosignals is the annotation of meaningful events that occur during the recording session, synchronously with the biosignal data. This ergonomic handheld manual button was specifically designed for that purpose.



Fig. 1. Ergonomic handheld design.




Fig. 2. Typical raw switch data (acquired with biosignals).

Pushbutton (BTN) Sensor Data Sheet

ORDERING GUIDE

Reference	Package Description
BTN1	Pushbutton (BTN) sensor in a standard configuration
BTN1-A-S	Pushbutton (BTN) sensor built with custom cable length A and custom sleeve color S; for standard physical characteristics in A or S use 0.